

The Axia SWC 72 Blanket Warmer

Features

Keypad temperature lock out

Push button digital operation

Large dual digital displays

**Displays actual temperature and set point
at all times**

**Separate digitally controlled heating
chambers**

**Rapid warm time (2-6 hours depending on
load)**

Uniform heating

Fully insulated

Audible and Visual Over-Temp Alarm

U.L./C.U.L. approved

FDA Registered

OSHDP Pre-Appro



Full Size Blanket Warmer

The **Axia SWC 72** full size blanket warmer is equipped with push button digital operation, large digital display, and separate digitally controlled heating chambers.

This full size blanket warmer is fully insulated and features uniform heating between both chambers. This warming cabinet operates rapidly with a warming time of 2 - 6 hours depending on load.



Axia SWC 72 Technical Specifications

Overall Size

Width	30 in (76.2 cm)
Depth	23 in (58.42 cm)
Height	74.5 in (189.23 cm)
Shelving	2 Fixed Shelves

Electrical

Power	120V, 7.5 amp., 50/60HZ .53 kw/hr, 1810 btu/hr UL/CUL listed FDA registered
-------	--

Chamber I.D.

Width	26 in (66.04 cm)
Depth	23 in (58.42 cm)
Height	61 in (154.94 cm)
Full Storage	21 Cubic Feet of Storage
Temperature Range	90-160 °F (32.2 - 71.1 °C)
Accuracy	± 3 °F of set point