

## Axia DWC 24 Blanket Warmer



### Features:

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L./C.U.L. approved
- FDA Registered
- OSHPD Pre-Appro

### Specifications

#### Overall Size

- 30"W x 26.5"D x 74.5"T
- 17.2 cubic feet of storage
- 3 adjustable shelves

#### Electrical

- 120V, 12 amp., 50/60HZ
- .75 kw/hr, 2557 btu/hr
- UL/CUL listed
- FDA registered

#### Upper Chamber I.D.

- 26"W x 23"D x 15.25"T
- 5.27 cubic feet of storage
- Temperature range: 90-160 F.
- Accuracy +/- 1 degree F. of set point
- 

#### Lower Chamber I.D.

- 26"W x 23"D x 34.5"T
- 11.93 cubic feet of storage
- Temperature range: 90-160 F.
- Accuracy +/- 3 degree F. of set point