

Axia Flex 1000R Surgical Table

Features:

Provides the features of low tabletop height adjustment for 520 to 1,010 (20.5" to 39.8", w/o pad) and tabletop 180 degree rotation for requirement of low height positioning needs in Neurosurgery or ENT, greater accessibility and upper body imaging without interference of the column base



Specifications

- **Elevation:** 20.5 ~ 39.8" (w/o mattress pad)
- **Head Section:** 60°/-90° (up/down)
- **Back Section:** 80°/40° (up/down)
- **Divided Leg Section:** 30°/90°, 0° ~ 90° (up/down, horizontal)
- **Trendelenburg/Reverse:** 26° each
- **Lateral Tilt:** Each 24° (left/right)
- **Flex / Reflex:** 220° / 100°
- **Table Top Rotation:** 360°
- **Auto-leveling function**