


Axia Flex 2000S Surgical Table

Features:

Provides the features of low tabletop longitudinal slide and extended height adjustment range from 685 to 1,135mm(27" - 44.7", w/o pad). The maximum tabletop height can fully meet the demanding specialties required by CRS surgery, Hip Replacement, L-Plating form intertrochanteric Fracture, Fixation of femoral head in orthopedics.



Specifications

- **Elevation:** 27" ~ 44.7"(w/o mattress pad)
- **Head Section:** 60°/-90° (up/down)
- **Back Section:** 80°/40° (up/down)
- **Divided Leg Section:** 30°/90°, 0° ~ 90° (up/down, horizontal)
- **Trendelenburg/Reverse:** 30° each
- **Lateral Tilt:** Each 20° (left/right)
- **Flex / Reflex:** 220° / 100°
- **Table Top Slide:** 7 1/2" ~ 13 3/4"
- **Table Top Rotation:** 180°
- **Auto-leveling function:** Returns all functions to level