

### Features:

- 530 mm of table top slide
- Comprehensive flexibility for patient positioning
- Position detecting alarm system
- Position memory system to facilitate surgery procedure



### Specifications

- **Table Top Dimensions:** 508mm W x 2000mm L ~ 2060 mm without side rails
  - **Elevation:** 680 ~1130 mm without mattress pad
  - **Head Section:** 60°/-90° up/down
  - **Back Section:** 80°/40° up/down
  - **Powered Leg Section:** 0°/100° up/down
  - **Trendelenburg/Reverse:** 25° / 30°
  - **Lateral Tilt:** Each 25° left/right
  - **Flex / Reflex:** 220° / 100°
  - **Table Top Slide:** 530 mm
  - **Powered Kidney Bridge Elevation:** 140mm
-