



Features:

- Excellent in Design and Performance
- Cutting Edge LED Technology (Truly White Light)
- Better Effect of Reflection with Unique Focusing Adjustment Function
- An Efficient Performance Completely Improve the Environments
- Integrated Surgical Lighting & Visualization System
- Flexible

Specifications:	Axia LED 500
Luminance (Lux)	80000*
Diameter of Light Head (CM)	50
Number of Bulbs	57
Color Temperature	4300
Color Render Index (Ra)	95
Color Render Index (R9)	94
Light Field Diameter (CM)	12-35
Focusing Depth (CM)	>90
Working Distance (CM)	50-80
Expected LED Service Life (H)	>30000

* The Luminance is optional to be customized from 80000 to 100000 Lux