



Features:

- Excellent in Design and Performance
- Cutting Edge LED Technology (Truly White Light)
- Better Effect of Reflection with Unique Focusing Adjustment Function
- An Efficient Performance Completely Improve the Environments
- Integrated Surgical Lighting & Visualization System
- Flexible

Specifications:	Axia LED 700	Axia LED 500
Luminance (Lux)	130000*	80000**
Diameter of Light Head (CM)	70	50
Number of Bulbs	108	57
Color Temperature	4300	4300
Color Render Index (Ra)	95	95
Color Render Index (R9)	94	94
Light Field Diameter (CM)	12-35	12-35
Focusing Depth (CM)	>90	>90
Working Distance (CM)	50-80	50-80
Expected LED Service Life (H)	>30000	>30000

* The Luminance is optional to be customized from 130000 to 160000 Lux

** The Luminance is optional to be customized from 80000 to 100000 Lux