


### Features:

The Axia STL4 is a mobile imaging table which is intended to support a 500 pound patient. The Axia STL4 is intended to be used with mobile x-ray for minimally invasive procedures.

### 4 Move Imaging Table:

- Motor Driven 29" to 39" Height Adjustment
- Motor Driven Lateral Tilt:  $\pm 17^\circ$
- Motor Driven Trendelenburg:  $\pm 13^\circ$
- Motor Driven Longitudinal Travel: 10"  $\pm 5$ "



### Specifications:

- Hi/Low
- Trendelenburg/Reverse Trendelenburg
- Lateral Tilt
- Longitudinal Travel
- 24" x 84" Carbon Fiber Diving Board
- 500 lb. Patient Weight Capacity
- 5" Locking Casters
- Hand Held Remote Control
- 120VAC, 60Hz., 6.5 Amp Power

### Dimensions:

- Table top: 84"L x 24"W x 1.62"T
  - Table dimensions: 86"L x 26"W
-