

Axia SWC 36 Blanket Warmer



Features:

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L./C.U.L. approved
- FDA Registered
- OSHPD Pre-Appro

Specifications

Overall Size

- 30"W x 26.5"D x 36"T
- 2 Adjustable Shelves

Chamber I.D.

- 26"W x 23"D x 25"T
- 8.63 cubic feet of storage
- Temperature range: 90-160 F.
- Accuracy +/- 2 F. degree of set point

Electrical

- 120V, 7.5 amp., 50/60HZ
- .35 kw/hr, 1195 btu/hr
- UL/CUL listed
- FDA registered