

Blanket Warmer



Axia DWC 24 Blanket Warmer

Displays

- Large dual digital displays
- Displays actual temperature and set point at all times
- Keyed temperature lockout

Heating

- Fully insulated
- Uniform heating
- Rapid Warm (2-6 hours)
- Separate digitally controlled heating chambers

Data Logging

- Easy Log integrated data logger with removable USB
- Automatically records data every 30 minutes
- Internal storage up to 2-years

Alarms

- Audible and visual Over-Temp alarms

Approximate Capacity

- 52 blankets & 120 solution bottles (1 Liter each)

Warranty

- 2-year warranty
-

Full Sized Warmer

Axia Blanket Warmers undergo stringent testing procedures during the manufacturing process to assure long life, the highest quality, and accurate temperature settings.

The Axia DWC series blanket warmers are the most user-friendly temperature recording devices on the market. They are equipped with independent, digitally controlled heating chambers that offer actual temperature and set point displays. Simply plug and play data requires no additional software.



Axia DWC 24 Technical Specifications

Overall Dimensions

Height	74.5" (189.2 cm)
Width	30" (76.2 cm)
Depth	26.5" (67.3 cm)
Cubic Feet Storage	17.21 cubic feet storage
Adjustable Shelves	3

Upper Chamber I.D.

Height	15.25" (38.73 cm)
Width	26" (66.04 cm)
Depth	23" (58.42 cm)
Cubic Feet Storage	5.28 cubic feet storage
Temperature Range	90°-160°F (32°-71°C)

Approximate Capacity

Aprx. Capacity	52 blankets & 120 solution bottles (1 Liter Each)
----------------	---

Electrical

Volts	120V
Amps	12.5 Amp
Hz	50/60 Hz
kWh (Avg)	.90 kWh
BTU/hr (Avg)	3071 BTU/hr

Lower Chamber I.D.

Height	34.5" (87.63 cm)
Width	26" (66.04 cm)
Depth	23" (58.42 cm)
Cubic Feet Storage	11.93 cubic feet storage
Temperature Range	90°-160°F (32°-71°C)