

### Axia UroMax 4 Surgical Table

#### Features

**Carbon fiber with integral facial cutout section**

**Leg extension easy to remove and very light weight**

**Diving board style for easy C-arm access**

**Large patient capacity (500-lb pound load)**

**A 2" patient high-quality comfort pad**


**Light weight – easy to physically move the table**



# Modern Surgical Table

Axia's **UroMax 4** C-arm surgical tables offer an excellent solution for real-time urological X-ray procedures. The cantilevered, carbon fiber, radiolucent top areas allow full fluoroscopic visualization and unobstructed C-arm positioning.

The UroMax 4 has highly versatile movements, options, and accessories. The carbon fiber leg extension accessory can quickly turn your Urology table into a procedure table for pain management and other C-arm guided procedures.



## Axia UroMax 3 Technical Specifications

### Dimensions

Table Length	55 in (139 cm)
Table Width	24 in (60.9 cm)
Table Height	30 - 40 in (76.2 - 101.6 cm)
Weight Capacity	500 lbs (226.7 kg)
Table Weight	740 lbs (335.7 kg)
H-Shaped Steel Base	51 x 26 in (129.5 - 66 cm)
Mattress	2 in thick (5 cm)
Radiolucent Area	27 x 24 in (68.5 x 60.9 cm)
Table Top	55 x 24 x 1.62 in (139 x 60 x 4 cm)

### Environmental

Storage Humidity	0- 92%
Operating Humidity	0-92%
Storage Temperature	32°F - 104°F (0°C - 40°C)
Operating Temperature	Up to 104°F (40°C)

### Tilts

Lateral Tilt	± 17°
Trendelenburg	± 13°
Longitudinal Travel	10 in (± 5 in)

### Power

AC Voltage	120V~
Frequency	60 Hz
AC Current	6.5A
Electric Shock Protection	Class I
Applied Part Rating	Type B
Mode of Operation	10% Duty Cycle 1 min on, 9 min off