

Axia SB1  
Imaging Table

Features

**Cantilevered (divingboard) style design**

**“H” style base**

**500 pound load rating**

**2” patient pad**

**24” wide top**

**Low boarding height**



# Modern Imaging Table

The **Axia SB1** imaging table is cantilevered (diving board) style design. The “H” style base allows maximum C-Arm clearance along the entire cantilevered length. It comes standard with 500 point load rating, 2” patient pad, 24” wide top, locking swivel castors, narrow base for easy patient transfer, low boarding height.



## Axia TL5 Technical Specifications

### Dimensions

Height	28 to 38” (71.12 to 96.52 cm)
Width	26” (66.04 cm)
Length	86” (218.44 cm)
Weight	520 lbs (236 kg)
Mattress Thickness	2” (5.08 cm)
Tabletop Length	84” (213.36 cm)
Tabletop Width	24” (60.96 cm)
Radiolucent Length	56” (142.24 cm)
Radiolucent Width	24” (60.96 cm)

### Base

Shape	H Shaped
Length	51” (129.54 cm)
Width	26” (66.04 cm)
Castors	5” with full swivel and wheel lock
Finish	Epoxy Finish

### Movements

Height	28 to 38” (71.12 to 96.52 cm)
--------	-------------------------------

### Power

AC Voltage	120V~
Frequency	60 Hz
AC Current	2.5A
Electric Shock Protection	Class I
Applied Part Rating	Type B
Mode of Operation	10% Duty Cycle 1 min on, 9 min off