


Hawk Series

Nothing Comes Close To this

**Axia
Surgical**



- Touchscreen intensity control
- Color temperature control
- Endoscopy mode
- Aluminum Dome
- Autoclavable handle
- Stainless steel or painted body

CEILING MOUNT & ROLLING STAND

Proven Technology & Trusted Quality

CE
2893

Axia Surgical
1.860.481.5383
166 Highland Park Drive
Bloomfield, CT 06002

axiasurgical.com
sales@axiasurgical.com

Hawk 917



Intensity: 160,000 LUX	LED Average Life >50,000 Hrs.
Field Diameter: 150-200 mm	Number of LEDs: 90
Depth Illumination 100-150 mm	Intensity Controller: Touchscreen with remote
Color Temp.: 3,800-4,800 K	Power Supply: 100-240V 50/60Hz AC
CR Index: 93 Ra	Power Consumption: 45W
Dome Diameter: 500 mm	

Hawk 918

Intensity: 160,000 LUX	LED Average Life >50,000 Hrs.
Field Diameter: 150-200 mm	Number of LEDs: 90
Depth Illumination 100-150 mm	Intensity Controller: Touchscreen with remote
Color Temp.: 3,800-4,800 K	Power Supply: 100-240V 50/60Hz AC
CR Index: 93 Ra	Power Consumption: 45W
Dome Diameter: 500 mm	



Hawk 919



Intensity: 160,000 LUX	LED Average Life >50,000 Hrs.
Field Diameter: 150-200 mm	Number of LEDs: 90 x 2
Depth Illumination 100-150 mm	Intensity Controller: Touchscreen with remote
Color Temp.: 3,800-4,800 K	Power Supply: 100-240V 50/60Hz AC
CR Index: 93 Ra	Power Consumption: 45W x 2
Dome Diameter: 500 mm	