

### Axia BariXray 5 Imaging Table

#### Features

**Cantilevered (diving board) style design**

**“H” style base**

**1000 pound load rating**

**2” patient pad**


**27” wide top**

**Low boarding height**



# Bariatric Imaging Table

The **Axia BariXray 5** imaging table is cantilevered (diving board) style design with bellows. The “H” style base allows maximum C-Arm clearance along the entire cantilevered length. It comes standard with 1000 pound load rating, 2” patient pad, 27” wide top, locking swivel castors, narrow base for easy patient transfer, extremely low boarding height.



## Axia BariXray 5 Technical Specifications

### Dimensions

Height	35 to 43” (88.9 to 109.2 cm)
Width	27” (68.6 cm)
Length	86” (218.4 cm)
Weight	875 lbs (396.9 kg)
Mattress Thickness	2” (5.08 cm)
Tabletop Length	84” (213.36 cm)
Tabletop Width	27” (68.6 cm)
Radiolucent Length	56” (142.2 cm)
Radiolucent Width	24” (61 cm)
Patient Capacity	1000 lbs (453 kg)

### Base

Shape	H Shaped
Length	63” (129.54 cm)
Width	26” (66.04 cm)
Castors	5” with full swivel and wheel lock
Finish	Epoxy Finish

### Movements

Height	35 to 43” (88.9 to 109.2 cm)
Lateral Tilt	±17°
Trendelenburg	±12°
Longitudinal Travel	16” (±5”) (40.6 cm (±12.7 cm)
Lateral Travel	10” (±5”) (25.4 cm (±12.7 cm)

### Power

AC Voltage	120V~
Frequency	60 Hz
AC Current	10A
Electric Shock Protection	Class I
Applied Part Rating	Type B
Mode of Operation	10% Duty Cycle 1 min on, 9 min off