

Axia STL5 Imaging Table

Features

2" patient high-quality comfort pad.

Facial cut-out for prone procedures.


Lightweight and easy to physically move the table.

24" table top which eliminates the need for arm boards.

500-pound patient weight capacity.

Hand pendant controls.

5 movement table



Modern Imaging Table

Axia's **STL5** is an X-ray imaging table is designed for high patient throughput for surgery centers, pain management, and other multidisciplinary medical facilities.

The Axia STL5 is a diving board style table with easy C-arm access. It has a large patient capacity of 500-pounds making it a great fit for management, minimally invasive spinal surgery, GI procedures, ultrasound procedures, orthopedic, endovascular, gall bladder, and other X-ray procedures.



Axia TL5 Technical Specifications

Dimensions

Height	33.5 to 43.5" (85.1 to 110.5 cm)
Width	26" (66 cm)
Length	86" (218.4 cm)
Weight	695 lbs (315 kg)
Mattress Thickness	2" (5.1 cm)
Tabletop Length	84" (213.4 cm)
Tabletop Width	24" (61 cm)
Radiolucent Length	56" (142.2 cm)
Radiolucent Width	24" (61 cm)

Base

Shape	I Shaped
Length	82" (208.3 cm)
Width	26" (66 cm)
Castors	5" with full swivel and wheel lock
Finish	Epoxy Finish

Movements

Height	33.5" to 43.5"
Lateral Tilt	±17°
Trendelenburg	±13°
Longitudinal Travel	±5"
Lateral Travel	10" (±5")

Power

AC Voltage	120V~
Frequency	60 Hz
AC Current	6.5A
Electric Shock Protection	Class I
Applied Part Rating	Type B
Mode of Operation	10% Duty Cycle 1 min on, 9 min off